Model 3900S

CW Ti:Sapphire Laser

User's Manual



The Solid-State Laser Company

1335 Terra Bella Avenue Mountain View, CA 94043

This manual contains information you need in order to safely install, align, operate, maintain, and service your *Model 3900S* CW Ti:sapphire laser.

The "Unpacking and Inspection" section contains information on how to unpack your Tsunami system and lists the items you should have received. It also lists any options you may have purchased. Please read this short, but important, section before you begin to unpack your unit.

The "Introduction" contains a brief description of the *Model 3900S* and lists system components and accessories available for this unit.

Following that section is an important chapter on laser safety. The *Model 3900S* and its pump laser are Class IV lasers and, as such, emit laser radiation which can permanently damage eyes and skin. This section contains information about these hazards and offers suggestions on how to safeguard against them. It also suggests installation procedures and maintenance you must perform in order to keep your system in compliance with CE and CDRH regulations. To ensure your system is installed according to these regulations and to minimize the risk of injury or expensive repairs, be sure to read this chapter—then follow these instructions.

"Laser Description" contains a brief exposition on Titanium:sapphire laser theory, which is followed by a more detailed description of the *Model 3900S* laser system. The chapter concludes with system specifications and outline drawings.

The middle chapters describe the *Model 3900S* controls, then guide you through its installation, alignment and operation. The last part of the manual covers maintenance and service. The latter includes a replacement parts list. "Customer Service" contains a list of world-wide Spectra-Physics service centers you can call if you need help.

Where-as the "Maintenance" section contains information you need to keep your laser clean and operational on a day-to-day basis, "Service and Repair" is intended to help you guide your Spectra-Physics field service engineer to the source of any problems. *Do not attempt repairs yourself while the unit is still under warranty*; instead, report all problems to Spectra-Physics for warranty repair.

Should you experience any problems with any equipment purchased from Spectra-Physics, or you are in need of technical information or support, please contact Spectra-Physics as described in "Customer Service." This chapter contains a list of world-wide Spectra-Physics Service Centers you can call if you need help.

This product has been tested and found to conform to "Directive 89/336/EEC for Electromagnetic Compatibility." Class A compliance was demonstrated for "EN 50081-2:1993 Emissions" and "EN 50082-1:1992 Immunity" as listed in the official *Journal of the European Communities*. It also meets the intent of "Directive 73/23/EEC for Low Voltage." Class A compliance was demonstrated for "EN 61010-1:1993 Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory use" and "EN 60825-1:1992 Radiation Safety for Laser Products." Refer to the "EC Declaration of Conformity" in Chapter 2, "Laser Safety."

Every effort has been made to ensure that the information in this manual is accurate. All information in this document is subject to change without notice. Spectra-Physics makes no representation or warranty, either express or implied, with respect to this document. In no event will Spectra-Physics be liable for any direct, indirect, special, incidental or consequential damages resulting from any defects in this documentation.

Finally, if you encounter any difficulty with the content or style of this manual, please let us know. The last page is a form to aid in bringing such problems to our attention.

Thank you for your purchase of Spectra-Physics instruments.

Environmental Specifications

CE Electrical Equipment Requirements

For information regarding the equipment needed to provide the electrical service listed under "Service Requirements" at the end of Chapter 3, please refer to specification EN-309, "Plug, Outlet and Socket Couplers for Industrial Uses," listed in the official *Journal of the European Communities*.

Environmental Specifications

The environmental conditions under which the laser system will function are listed below:

Indoor use

Altitude: up to 2000 m Temperatures: 10° C to 40° C

Maximum relative humidity: 80% non-condensing for temperatures up to

31° C.

Mains supply voltage: do not exceed $\pm 10\%$ of the nominal voltage

Insulation category: II Pollution degree: 2

Table of Contents

Preface	iii
Environmental Specifications	
CE Electrical Equipment Requirements	
Warning Conventions	xi
Standard Units	xiii
Unpacking and Inspection	xv
Unpacking Your Laser	
Chapter 1: Introduction	1-1
System Components	
Chapter 2: Laser Safety	2-1
Precautions for the Safe Operation of Class IV-High Power Lasers. Interlocks Maintenance Required to Keep this Laser Product in Compliance Radio logical Health (CDRH) Regulations CE/CDRH Radiation Control Drawing Label Translations CE Declaration of Conformity (emissions) CE Declaration of Conformity (low voltage) Sources for Additional Information Laser Safety Standards Equipment and Training	2-3 re with Center for Devices and
Chapter 3: Laser Description	3-1
Ti:sapphire as a Laser Medium Four-Mirror Folded Cavity Astigmatic Correction and Pumping Birefringent Filter Tuning Characteristics of the CW Ti:sapphire Laser Optics Sets Pumping Options	3-2 3-4 3-4 3-5

Chapter 4: Installation and Alignment	
Controls and Connections	
Cavity Alignment	2
Chapter 5: Narrow Linewidth Options5-1	
Thin Etalon Assembly 5-1 Installation 5-1 Alignment 5-2 Thick Etalon Assembly 5-5 Installation 5-5 Alignment 5-6	
Chapter 6: Operation6-1	
System Startup	
Chapter 7: Maintenance	
Notes on the Cleaning of Laser Optics 7-1 Equipment Required 7-2 Cleaning Solutions Required 7-3 General Procedures for Cleaning Optics 7-3 Steps for Cleaning Optics 7-5 Removing the Mirrors for Cleaning 7-5 Output Mirror M7 7-5 High Reflector M8 7-5 Fold Mirrors M5 and M6 7-5 Pump Mirrors M1, M2, M3 and M4 7-6 Cleaning the Ti:sapphire Laser Rod Surfaces 7-6 Cleaning the Birefringent Filter 7-6 Cleaning the Brewster Windows 7-7 Replacing Filters in the Optional Model 3910 Purge Unit 7-7	
Chapter 8: Service and Repair	
Troubleshooting Guide	
Chapter 9: Customer Service9-1	
Warranty	

Appendix A: Material Safety Data Sheets (MSDS)	l
Notes	
Report Form for Problems and Solutions	
List of Figures	
Figure 1-1: The Model 3900S Ti:sapphire Laser Figure 2-1: These CDRH and CE standard safety warning labels would be appropriate for use as entry warning signs (ANSI 4.3.10.1, EN 60825-1). Figure 2-2: Folded Metal Beam Target 2-2: Figure 2-3: The Model 3900S Interlock Shutter 2-3: Figure 2-4: Model 3900S CE/CDRH Radiation Control Drawing 2-5: Figure 3-1: Energy level structure of Ti3+ in sapphire 3-1: Energy level structure of Ti3+ in sapphire 3-1: Figure 3-2: Absorption and emission spectra of Ti:sapphire 3-2: Figure 3-3: Folded four-mirror cavity with longitudinal pumping 3-3: Figure 3-4: Typical three-mirror cavity used in a CW dye laser. 3-3: Figure 3-4: Typical three-mirror cavity used in a CW dye laser. 3-3: Figure 4-1: Rod Translation Stage Lock-down Screws 4-2: Model 3900S Optical Element Identification 4-2: Figure 4-2: Model 3900S Optical Element Identification 4-2: Figure 4-3: Model 3900S Mirror Mount and Control Locations 4-3: Figure 4-4: Model 3910 Regulator/Filter Purge Unit 4-4: Figure 4-5: Typical system layout for Model 3900S and pump laser. 4-5: Figure 4-6: Location of M5 and M6 lock-down screws 4-8: Figure 4-7: M4 Lock-down Screws and Micrometer Setting 4-8: Figure 4-8: The fluorescence spot as it appears on mirror M7. 4-11: Figure 4-9: Wavelength Calibration Chart 4-15: Figure 4-9: Wavelength Calibration Chart 4-16: Figure 5-1: Model 3900S Mirror Mount and Control Locations, including Etalon Assembly Options 5-3: Figure 5-2: Model 3900S Mirror Mount and Control Locations, including Etalon Assembly Options 5-3: Figure 5-3: Etalon Mount Assembly Options 5-3: Figure 5-3: Etalon Mount Assembly Options 5-3: Figure 5-3: Etalon Mount Assembly With M8 Modification 5-4 Figure 7-2: Lens Tissue Folded for Cleaning 7-4 Figure 7-3: Underside of the Model 3910, showing filter assembly placement and orientation 7-8	0 3 3 3
List of Tables	
Table 2-1 : Label Translations.2-6Table 3-1 : Model 3900S Specifications3-6Table 8-1 : Replacement Parts, Standard Model 3900S Optics8-2Table 8-2 : Model 3900S Optic Sets8-2Table 8-3 : Model 3900S Accessories and Replacement Parts8-3Table 8-4 : Hardware8-3	

Warning Conventions

The following warnings are used throughout this manual to draw your attention to situations or procedures that require extra attention. They warn of hazards to your health, damage to equipment, sensitive procedures, and exceptional circumstances. All messages are set apart by a thin line above and below the text as shown here.



Laser radiation is present.

Danger!



Condition or action may present a hazard to personal safety.





Condition or action may present an electrical hazard to personal safety.

Warning!



Condition or action may cause damage to equipment.

Warning! ESD



Action may cause electrostatic discharge and cause damage to equipment.





Condition or action may cause poor performance or error.

Note



Text describes exceptional circumstances or makes a special reference.

Don't Touch!



Do not touch.

Eyewear Required



Appropriate laser safety eyewear should be worn during this operation.



Refer to the enclosed documents and manual before operating or using this device.

Standard Units

The following units, abbreviations, and prefixes are used in this Spectra-Physics manual:

Quantity	Unit	Abbreviation
mass	kilogram	kg
length	meter	m
time	second	S
frequency	hertz	Hz
force	newton	N
energy	joule	J
power	watt	W
electric current	ampere	Α
electric charge	coulomb	С
electric potential	volt	V
resistance	ohm	Ω
inductance	henry	Н
magnetic flux	weber	Wb
magnetic flux density	tesla	Т
luminous intensity	candela	cd
temperature	celcius	С
pressure	pascal	Pa
capacitance	farad	F
angle	radian	rad

			P	refixes				
tera	(1012)	Т	deci	(10-1)	d	nano	(10-9)	n
giga	(10°)	G	centi	(10-2)	С	pico	(10-12)	р
mega	(10^6)	М	mill	(10-3)	m	femto	(10-15)	f
kilo	(10^3)	k	micro	(10-6)	μ	atto	(10-18)	а

Unpacking and Inspection

Unpacking Your Laser

Your laser was packed with great care and its containers inspected prior to shipment. It left Spectra-Physics in good condition. Upon receiving your laser, immediately inspect the outside of the shipping containers. If there is any major damage (holes in the containers or cracked wooden frame members), insist that a representative of the carrier be present when you unpack the contents.

Carefully inspect your laser as you unpack it. If you notice any damage, such as dents or scratches on the cover, or broken knobs, immediately notify the carrier and your Spectra-Physics sales representative.

Keep the shipping containers. If you file a damage claim, you may need the containers to demonstrate that the damage occurred as a result of shipping. If you need to return the laser for service, these specially designed crates assure adequate protection.

Warning!



Spectra-Physics considers itself responsible for the safety, reliability, and performance of the *Model 3900S* only under the following conditions:

- All field installable options, modifications, or repairs are performed by persons trained and authorized by Spectra-Physics.
- The equipment is used only according to instructions provided in this manual.

System Components

- Model 3900S laser head
- *Model 3910* nitrogen purge unit (optional)

Accessory Kit

Included with the laser system is this manual, a test summary, a packing slip listing all the components shipped with this order, and an accessory kit containing the following items:

- Two pump beam routing mirror assemblies and a set of dust tubes
- Two hoses (cooling system water supply and return lines)
- Three foot clamps with hardware for mounting the *Tsunami* laser head

- A tool kit containing:
 - Allen (hex) ball drivers for optimizing laser output
 - An IR card
 - a tweezer
 - A plastic hemostat
 - A packet of Kodak Lens Cleaning Paper[™]
 - Any optional optics and optical assemblies that were ordered

You will need to supply several items, including:

- Spectrophotometric-grade (HPLC) acetone and methanol for optics cleaning
- Clean, lint-free finger cots or powderless latex gloves for optics cleaning
- Ir viewer
- Optical spectrum analyzer (recommended)
- Fast photodiode (recommended)

Chapter 1 Introduction

Titanium-doped sapphire (Ti:sapphire) is a solid-state laser medium capable of tunable operation over a broad range of near infrared (IR) wavelengths.

Because of Ti:sapphire's broad absorption band in the blue and green, energy for the lasing process (the "pump" energy) can be supplied by a *Millennia*® series continuous wave (CW), diode-pumped, solid-state laser or by a *Model 2016, 2017, 2020, 2030, 2040, 2060* and *2080* argon ion laser.

Because there are no dyes to change, no circulators to clean, and no drop in output power as dyes degrade, Ti:sapphire solid-state lasers eliminate much of the effort normally required to operate dye lasers in the near infrared.

System Components

- Model 3900S laser head
- *Model 3910* nitrogen purge unit (optional)



Figure 1-1: The Model 3900S Ti:sapphire Laser

Patents

Danger!



The Spectra-Physics *Model 3900S* and its pump laser are Class IV-High Power Lasers, whose beams are, by definition, safety and fire hazards. Take precautions to prevent accidental exposure to both direct and reflected beams. Diffuse as well as specular reflections can cause severe eye or skin damage.



Because the *Model 3900S* laser emits cw infrared radiation, it is extremely dangerous to the eye. Infrared radiation passes easily through the cornea, which focuses it on the retina, where it can cause instantaneous permanent damage.

Precautions for the Safe Operation of Class IV-High Power Lasers.

Eyewear Required



- Wear protective eyewear at all times; selection depends on the wavelength and intensity of the radiation, the conditions of use, and the visual function required. Protective eyewear vendors are listed in the Laser Focus World, Lasers and Optronics, and Photonics Spectra buyer's guides. Consult the ANSI or ACGIH standards listed at the end of this section for guidance.
- Maintain a high ambient light level in the laser operation area. This
 keeps the eye's pupil constricted, thus reducing the possibility of eye
 damage.
- Keep the protective cover on the laser at all times.
- Avoid looking at the output beam; even diffuse reflections are hazard-
- Avoid wearing jewelry or other objects that may reflect or scatter the beam while using the laser.
- Use an infrared detector or energy detector to verify that the laser beam is off before working in front of the laser.
- Operate the laser at the lowest beam intensity possible, given the requirements of the application.
- Expand the beam whenever possible to reduce beam power density.
- Avoid blocking the output beam or its reflection with any part of your body.

- Establish a controlled access area for laser operation. Limit access to those trained in the principles of laser safety.
- Post prominent warning signs near the laser operation area (Figure 2-1).
- Set up experiments so the laser beam is either above or below eye level.
- Provide enclosures for beam paths whenever possible.
- Set up shields to prevent specular reflections.
- Set up an energy absorbing target to capture the laser beam, preventing unnecessary reflections or scattering (Figure 2-2).



Figure 2-1: These CDRH and CE standard safety warning labels would be appropriate for use as entry warning signs (ANSI 4.3.10.1, EN 60825-1).




Figure 2-2: Folded Metal Beam Target





Use of controls or adjustments, or performance of procedures other than those specified herein may result in hazardous radiation exposure.

Follow the instructions contained in this manual for safe operation of your laser. At all times during operation, maintenance, or service of your laser, avoid unnecessary exposure to laser or collateral radiation* that exceeds the accessible emission limits listed in "Performance Standards for Laser Products," *United States Code of Federal Regulations*, 21CFR1040 10(d).

^{*} Any electronic product radiation, except laser radiation, emitted by a laser product as a result of, or necessary for, the operation of a laser incorporated into that product.

Interlocks

Because the energy to drive the lasing process in the *Model 3900S* comes from another laser and not from an internal source (such as electrical discharge), the interlock differs slightly from that of other lasers.

The *Model 3900S* laser has only one interlock: the shutter interlock (Figure 2-3). It blocks the pump beam at the entrance to the laser head housing to prevent the *Model 3900S* from lasing. When installed, the laser head cover holds the shutter interlock open for normal operation. When the cover is removed, the shutter closes automatically.



Operating the laser with the cover off may expose people to high voltages and high levels of radiation. It also increases the rate of optical surface contamination. Therefore, operating the laser with the cover off is not recommended.

The alignment procedures in this manual require internal adjustments while the laser is operating. The interlock can be defeated to allow this. When the cover is removed and access to the pump beam is required, raise the red shutter lever to defeat the interlock and hold the shutter open. In this position, the red lever clearly shows the defeat and prevents cover installation until the shutter lever is lowered to the closed position.





Figure 2-3: The Model 3900S Interlock Shutter

Maintenance Required to Keep this Laser Product in Compliance with Center for Devices and Radio logical Health (CDRH) Regulations

This laser product complies with Title 21 of the *United States Code of Federal Regulations*, chapter 1, subchapter J, parts 1040.10 and 1040.11, as applicable. To maintain compliance with these regulations, once a year, or whenever the product has been subjected to adverse environmental conditions (e.g., fire, flood, mechanical shock, spilled solvent, etc.), check to see that all features of the product identified below function properly. Also, make sure that all warning labels remain firmly attached (refer to the CDRH/CE drawing later in this chapter).

- 1. Verify that removing the laser head cover closes the laser shutter and blocks the pump beam, thereby turning off the Ti:sapphire laser.
- 2. Verify that, when the cover interlock is defeated, the defeat mechanism is clearly visible and prevents installation of the cover until disengaged.
- 3. Verify that all labels listed in Figure 2-4, "Model 3900S Radiation Control Drawing," are present and are firmly affixed.
- 4. Verify removing the remote interlock plug on the pump laser prevents laser operation.
- 5. Verify the laser system will only operate when the pump laser's interlock key switch is in the ON position, and that the key can only be removed when the switch is in the OFF position.
- 6. Verify the emission indicator on the pump laser works properly; that is, it emits a visible signal whenever the laser is on.
- 7. Verify that the time delay between turn-on of the pump laser emission indicator and that starting of that laser gives you enough warning to allow action to avoid exposure to laser radiation.
- 8. Verify removing the cover of the pump laser shuts off the laser.
- 9. Verify, when the cover interlock on the pump laser is defeated, the defeat mechanism is clearly visible and prevents installation of the cover until disengaged.

CE/CDRH Radiation Control Drawing






Figure 2-4: Model 3900S CE/CDRH Radiation Control Drawing

Label Translations

For safety, the following translations are provided for non-English speaking personnel. The number in parenthesis in the first column corresponds to the label number listed on the previous page.

Table 2-1: Label Translations

1 -1 1 "	F '	0	0	D
Label #	French	German	Spanish	Dutch
(2) Caution, Defeat- able Interlock (EMI)	Attention. Rayonnement visible et invisible dangereux en cas d'ouverture et lorsque la sécurité est neutralisée.	Achtung! Sichtbare und unsichtbare schädliche elektromagnetische Strahlung wenn Abdeckung geöffnet und Sicherheitsverriegelung überbrückt. Bedienungsanleitung beachten!	Precaución, radiación peligrosa electromagnética visible e invisible con el dispositivo de seguridad abierto o con su indicación alterada.	Let op. Zichtbare en onzichtbare gevaar- lijke electromagne- tische straling indien geopend en interlock overbrugd.
(3) Danger, Defeat- able Interlock	Attention- Rayonnement Laser visible et invisible en cas D'Ouverture et lorsque la securite est neutralisse; exposition dangereuse de l'oeil ou de la peau au rayonnement dirct ou diffus.	Vorsicht; Austritt von sichtbarer un unsichtbarer un unsichtbarer Laserstruhlung, wenn Abdeckung geoffnet und Sicherhetisschalter uberbruckt; Bestrahlung von Auge oder Haut durch direkte oder Streustreustrahlung vermeiden.	Peligro, al abrir y retirer el dispositivo de seguridad exist radiacion laser visible e invisible; evite que los ohos o la piel queden expuestos tanto a la radiacion direta como a la dispersa.	Gevaar; zichtbare en niet zichtbare laser-straling wanneer geopend en bij uitgeschakelde interlock; Vermijd blootstelling van oog of huid aan directe straling of weerkaatsingen daarvan.
(4) Aperture Label	Ouverture Laser - Exposition Dan- gereuse - Un Rayon- nement laser visible et invisible est emis par cette ouverture.	Austritt von sichtbarer und unsictbarer Laserstrahlung; nicht dem Strahl aussetzen.	Por esta abertura se emite radiacion laser visible e invisible; evite la exposicion.	Vanuit dit apertuur wordt zichtbare en niet zichtbare laser- straling geemiteerd; vermijd blootstellilng.
(5) CDRH Logo- type Danger Label	Attention— Rayonnement Laser Visible et Invisible en Cas D'Ouverture et lorsque la securite est neutralisee; expo- sition dangereuse de l'oeil ou de la peau au rayonnement direct ou diffus. Puissance et longueurs D'onde dependant de la con- figuration et de la puissance de pompe. Laser de Classe 4.	Vorsicht; Austritt von sichtbarer un unsichtbarer un unsichtbarer Laserstrahlung wenn Abdeckung geoffnet und Sicherhetisschalter uberbruckt; Bestrahlung von Auge oder Haute durch direkte oder Streustrahlung vermeiden. Leistung, Wellenlange und Pulsbreite sind abhangig von Pumpquelle und Laserkonfiguration. Laserklasse 4.	Peligro, al abrir y retiar el dispositivo de seguridad exist radiacion laser visible e invisible; evite que los ohos o la piel queden expuestos tanto a la radiacion directa como a la dispersa. Potencia, Longitud de onda y anchura de pulso dependen de las opciones de bombeo y de la configuracion del laser. Producto laser clase 4.	Gevarr, zichtbare en neit zichtbare laserstraling wanneer geopend en bij uitgeschakelde interlock; Vermijd blootstelling van oog of huid aan directe straling of weerkaatsingen daarvan. Vermogen golfleugten en pulsduur afhankelijk van pomp optics en laser configuratie. Klasse 4 Laser Produkt.

CE Declaration of Conformity

We,

Spectra-Physics, Inc. Industrial and Scientific Lasers 1330 Terra Bella Avenue P.O. Box 7013 Mountain View, CA. 94039-7013 United States of America

declare under sole responsibility that the:

Model 3900S cw Ti:sapphire Laser,

Manufactured after December 31, 1996

meets the intent of "Directive 89/336/EEC for Electromagnetic Compatibility."

Compliance was demonstrated (Class A) to the following specifications as listed in the official *Journal of the European Communities*:

EN 50081-2:1993 Safety Requirements for Emissions:

EN 55011 Class A Radiated EN 55011 Class A Conducted

EN 50082-1:1992 Immunity:

IEC 801-2 Electrostatic Discharge IEC 801-3 RF Radiated IEC 801-4 Fast Transients

I, the undersigned, hereby declare that the equipment specified above conforms to the above Directives and Standards.

Steve Sheng

Vice President and General Manager

Spectra-Physics, Inc.

New Shong

Industrial and Scientific Lasers

January 1, 1998

CE Declaration of Conformity

We,

Spectra-Physics, Inc. Industrial and Scientific Lasers 1330 Terra Bella Avenue P.O. Box 7013 Mountain View, CA. 94039-7013 United States of America

declare under sole responsibility that the

Model 3900S cw Ti:Sapphire Laser,

Manufactured after December 31, 1996

meets the intent of "Directive 73/23/EEC, the Low Voltage directive."

Compliance was demonstrated to the following specifications as listed in the official *Journal of the European Communities*:

EN 61010-1: 1993 Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory use:

EN 60825-1: 1994 Safety for Laser Products.

I, the undersigned, hereby declare that the equipment specified above conforms to the above Directives and Standards.

Steve Sheng

Vice President and General Manager

Spectra-Physics, Inc.

Seno Shong

Industrial and Scientific Lasers

January 1, 1998

Sources for Additional Information

The following are some sources for additional information on laser safety standards, safety equipment, and training.

Laser Safety Standards

Safe Use of Lasers (Z136.1: 1993)

American National Standards Institute (ANSI)

11 West 42nd Street

New York, NY 10036

Tel: (212) 642-4900

Occupational Safety and Health Administration (Publication 8.1-7)

U. S. Department of Labor

200 Constitution Avenue N. W., Room N3647

Washington, DC 20210

Tel: (202) 693-1999

A Guide for Control of Laser Hazards, 4th Edition, Publication #0165

American Conference of Governmental and

Industrial Hygienists (ACGIH)

1330 Kemper Meadow Drive

Cincinnati, OH 45240

Tel: (513) 742-2020

Internet: www.acgih.org/home.htm

Laser Institute of America

13501 Ingenuity Drive, Suite 128

Orlando, FL 32826

Tel: (800) 345-2737

Internet: www.laserinstitute.org

Compliance Engineering

70 Codman Hill Road

Boxborough, MA 01719

Tel: (978) 635-8580

International Electrotechnical Commission

Journal of the European Communities

EN60825-1 TR3 Ed.1.0—Laser Safety Measurement and Instrumentation

IEC-309—Plug, Outlet and Socket Coupler for Industrial Uses

Tel: +41 22-919-0211

Fax: +41 22-919-0300

Internet: http://ftp.iec.c.h/

Cenelec

European Committee for Electrotechnical Standardization

Central Secretariat

rue de Stassart 35

B-1050 Brussels

Document Center

1504 Industrial Way, Unit 9

Belmont, CA 94002-4044

Tel: (415) 591-7600

Equipment and Training

Laser Safety Guide Laser Institute of America 12424 Research Parkway, Suite 125 Orlando, FL 32826 Tel: (407) 380-1553

Laser Focus World Buyer's Guide Laser Focus World Penwell Publishing 10 Tara Blvd., 5th Floor Nashua, NH 03062 Tel: (603) 891-0123

Lasers and Optronics Buyer's Guide Lasers and Optronics Gordon Publications 301 Gibraltar Drive P.O. Box 650 Morris Plains, NJ 07950-0650 Tel: (973) 292-5100

Photonics Spectra Buyer's Guide Photonics Spectra Laurin Publications Berkshire Common PO Box 4949 Pittsfield, MA 01202-4949 Tel: (413) 499-0514 Titanium doped sapphire (Ti:sapphire) is a solid-state laser medium capable of tunable laser action over a broad range of near infrared (IR) wave lengths. The broad absorption band in the blue and green make a *Millennia*® diode-pumped laser an ideal pump source for this material. Ti:sapphire eliminates much of the effort required to operate dye lasers in the near infrared, as there are no dyes to change, no circulators to clean, and no drop in output power as dyes degrade.

Ti:sapphire as a Laser Medium

Ti:sapphire is a crystalline material produced by introducing Ti₂O₃ into a melt of Al₂O₃. A boule of single crystal material is grown from this melt by one of several techniques. The Ti³⁺ion is responsible for the laser action of Ti:sapphire. It substitutes for a small percentage of the Al³⁺ ion, and the electronic ground state of the Ti³⁺ ion is split into a pair of vibrationally-broadened levels as shown in Figure 3-1.




Figure 3-1: Energy level structure of Ti³⁺ in sapphire

Absorption transitions occur over a broad range of wavelengths from 400 to 600 nm, only one of which is shown in Figure 3-1. Fluorescence transitions occur from the lower vibrational levels of the excited state to the upper vibrational levels of the ground state. The resulting emission and absorption spectra are shown in Figure 3-2.




Figure 3-2: Absorption and emission spectra of Ti:sapphire

The fluorescent lifetime of the excited states decrease due to non-radiative transitions as the temperature rises, therefore the material must be maintained at or below room temperature for efficient operation.

Although the fluorescence band extends from wavelengths as short as 600 nm to wavelengths greater than 1 mm, the long wavelength side of the absorption band overlaps the short wavelength end of the fluorescence spectrum. Therefore, laser action is only possible at wavelengths longer than 650 nm. The tuning range is further reduced by an additional weak absorption band that overlaps the fluorescence spectrum. This band has been traced to the presence of Ti⁴⁺ ions, but it is also dependent on material growth techniques and Ti³⁺ concentration. Finally, the tuning range is also affected by mirror coatings, tuning element losses, available pump power, and pump mode quality.

Four-Mirror Folded Cavity

For continuous-wave (CW) pumping, there is one basic requirement for laser action: the unsaturated round-trip cw gain must exceed the round-trip loss from all sources. The CW gain is obtained by having a high inversion density and a sufficient length of Ti:sapphire material. A high inversion density comes from having a high pump intensity and a high Ti³⁺ ion concentration. Losses in the Ti:sapphire laser come from losses in mirror coatings and polished surfaces and more importantly, in the residual loss in the Ti:sapphire material itself. This loss is both proportional to the rod length and varies with the Ti³⁺ concentration, generally increasing as the Ti³⁺ concentration increases.

It is very difficult to achieve a continuous, high inversion density over the entire volume of a Ti:sapphire rod several millimeters in diameter. To cir-

cumvent this problem, the pump light is focused to a narrow line within the rod, and the oscillating laser mode is similarly focused and overlapped within the same volume—a technique known as longitudinal pumping. A proprietary computer analysis developed at Spectra-Physics was used to optimize the optical cavity design of the Ti:sapphire laser. The result is a laser that is efficient, has a low pumping threshold, and is also easy to align.

A four-mirror folded cavity is used to allow tightly focused cavity and pump modes within the laser rod, while providing a collimated output beam of normal size. The beam path is shown in Figure 3-3.




Figure 3-3: Folded four-mirror cavity with longitudinal pumping.




Figure 3-4: Typical three-mirror cavity used in a CW dye laser.

The pump illumination in the Ti:sapphire laser must be collinear to the cavity mode over a relatively long length of laser rod. If the typical three-mirror cavity used in dye lasers (Figure 3-4) were used with Ti:sapphire, the pump and cavity modes would overlap at one of the fold mirrors and would make laser alignment difficult. In the *Model 3900S*, an additional collimated leg is used and the remaining pump beam is "dumped" through one of the fold mirrors. The resulting cavity is easy to align.

Astigmatic Correction and Pumping

Using a focusing mirror at an angle other than normal incidence leads to astigmatism unless corrected by some other element, e.g. the Brewster-angle Ti:sapphire rod. Folded cavities in which this astigmatism is not eliminated yield oval-shaped output beams that are difficult to focus. By a careful choice of the cavity fold angles and rod length, astigmatism in the output beam can be completely eliminated. But because astigmatism still exists within the laser rod, the pump beam must also be astigmatic to allow efficient coupling between the pump and intra-cavity beam. A concave focusing mirror used away from normal incidence at an angle optimized by a patented computer analysis allows astigmatism to be induced in the pump beam that most efficiently matches the astigmatism of the cavity mode. The resulting laser has high conversion efficiency and can be pumped with a variety of solid-state and argon ion pump lasers.

Birefringent Filter

The birefringent filter for the *Model 3900S* consists of three crystalline quartz plates placed within the laser cavity at Brewster's angle. The plates are cut parallel to their optical axes, and their birefringence causes the linear polarization of the incident laser beam to become elliptical. Only one frequency will make a complete 180° (or multiple there of) polarization flip and return to linear polarization; the polarization of all other beams rotates more or less than that and remains elliptically polarized. The elliptically polarized wavelengths suffer additional losses at all of the Brewsterangle surfaces within the cavity and fail to reach lasing threshold.

The birefringent filter used in the *Model 3900S* has a free spectral range a little greater than 150 nm. This free spectral range is the difference between the adjacent eigenwavelengths—those wavelengths that return to linear polarization after traversing the filter. Rotating the filter about an axis normal to the plates changes these eigenwavelengths.

A single birefringent filter is used for all the wavelengths the *Model 3900S* can access. Although the broad tuning range of Ti:sapphire over laps two "orders" of the filter, the narrower wavelength range of the mirror set permits only one order to oscillate at a time. The filter is installed and adjusted for the 700 to 850 nm range. If the 850 to 1000 nm range is desired, the only change required of the filter (along with a change in the optics set) is to rotate the entire three-plate filter in its holder to access the order most useful for longer wavelengths. The alignment of the three-plate stack itself is set at the factory and will never need changing.

Tuning Characteristics of the CW Ti:sapphire Laser

The tuning characteristics of the CW Ti:sapphire laser are influenced by several factors, some dependent on cavity parameters such as birefringent tuning elements and Brewster surfaces, and some inherent to the Ti:sapphire material itself.

Because both the Ti:sapphire and the filter plates are birefringent, proper filter operation is achieved only if the c-axis of the rod is aligned coplanar with the polarization of the electric field within the cavity. Since the birefringent filter plates and the Ti:sapphire rod represent a total of eight Brewster-angle surfaces, the polarization within the cavity is largely determined by the orientation of these surfaces. Cavity losses are minimized and tuning is optimized if all these surfaces are accurately aligned at Brewster's angle. The *Model 3900S* uses a proprietary Ti:sapphire rod holder (patent pending) to both orient the rod surfaces at Brewster's angle and to align the c-axis to be coplanar to the electric field vector. This technique is used to completely compensate for unavoidable errors in rod orientation during cutting and polishing.

Secondly, since the tuning range of Ti:sapphire is far broader than any other comparable infrared laser medium, the quartz plates used in the bire-fringent filter must be fabricated to much higher precision than is typically required for other tunable media such as infrared dyes.

Finally, the presence of airborne water vapor in the laser cavity leads to discontinuities in the tuning range of Ti:sapphire. Purging the unit with dry nitrogen gas eliminates most of the water present to improve the overall tuning performance of the system.

Optics Sets

The *Model 3900S* is supplied with two optics sets. The first allows tuning between 700 and 850 nm: the second between 850 and 1000 nm. Only the four cavity mirrors will require replacement when switching between these two tuning ranges; the pump mirrors remain unchanged for the entire tuning range of the *Model 3900S*. The mirrors are quickly and easily changed. Two of the optional optics sets available provide tunability from 675 to 750 nm and from 950 to 1100 nm. Thus, the total possible tuning range available from the *Model 3900S* with all four optics sets and a 20 W pump source is 675 to 1100 nm. The complete tuning range may not be available if less than 20 W of input pump power is used.

A third optional mirror set covers the 800 to 900 nm range. A special birefringent filter is required for use with this mirror set to avoid the flip between tuning orders with the standard birefringent filter.

Pumping Options

The *Model 3900S* can be pumped by Spectra-Physics *Millennia Vs J, VIs J, VIIIs J* sand Xs J solid-state lasers or *Model 2016, 2017, 2020, 2030, 2040, 2060* and 2080 argon ion lasers. The ion pump lasers can be operated single line at 514.5 nm or all-lines, as long as at least 5 W of TEM₀₀ output power is obtained. The most efficient operation is obtained with pump powers over 10 W. Call your Spectra-Physics service engineer regarding pumping options.

Table 3-1: Model 3900S Specifications¹

Performance ²		
Category	3900S-1	3900S-2
Average Power 5 W Pump 20 W Pump	750 mW 2.5 W	750 mW 3.25 W
Tuning Range 5 W Pump 20 W Pump	750–950 nm 700–980 nm	750–950 nm 700–980 nm
Peak Wavelength	790 nm	790 nm
Linewidth	<40 GHz	<40 GHz
Noise ³	<1%	<1%
Stability ⁴	<3%	<3%
Spatial Mode	TEM ₀₀	TEM ₀₀
Polarization	>100:1 Horizontal	>100:1 Horizontal
Beam Diameter ^{2,5}	0.95 mm	0.95 mm
Beam Divergence ^{2,6}	1.0 mrad	1.0 mrad

Specifications subject to change without notice. They are valid only if the pump laser meets all of its published specifications, is in POWER or LIGHT mode, and pump power does not exceed the minimum required to meet specified output

power.

² Specifications apply to operation at the peak wavelength noted. Contact your field sales engineer for specifications at specifications apply to operation at the peak wavelength noted. Contact your field sales engineer featened wavelengths.

RMS, measured in a 10 Hz to 1 MHz bandwidth.

Percent drift in any 2-hour period after a 1-hour warm-up and less than ±3° C temperature change.

At \frac{1}{2} points.

Full angle.

Mechanical	
Water	
Pressure Range	
Minimum	103 kPa (15 psi)
Maximum	517 kPa (75 psi)
Inlet Temperature Range	
Minimum	18°C (65°F)
Maximum	27°C (80°F)
Physical Dimensions	
Size (L x W x H)	65.1 x 24.1 x 20.6 cm (25.6 x 9.50 x 8.1 in.)
Weight	14.1 kg (31.0 lb)

Installation and Alignment

Warning!



The following installation and alignment procedures are to be performed by an authorized Spectra-Physics service engineer. Please call your service representative to arrange an installation appointment, which is part of your purchase agreement. Allow only those qualified and authorized by Spectra-Physics to install and set up your laser system. You will be charged for any damage incurred if you attempt the initial installation yourself, and such action may void your warranty. These instructions are provided in the event that you move the laser to a new location in the future.

Controls and Connections

Model 3900S Laser Head

The faces of the mirror control knobs have been color-coded for your convenience: vertical is blue, horizontal is green.

Polarization rotator mirrors (M₁ and M₂)—rotate the vertically polarized pump beam to horizontal polarization and turn the beam by 90°.

Pump beam mirror (M_3) —directs the pump beam onto the center of the pump focus mirror M_4 .

Pump focus mirror (M₄)—directs and focuses the pump beam into the Ti:sapphire rod. Both its orientation and distance from the Ti:sapphire rod are adjustable.

Cavity fold mirrors (M₅ and M₆)—focus the cavity mode in the Ti:sapphire rod while transmitting the pump beam.

Output mirror (M_7) and cavity high reflector (M_8) —complete the four-mirror cavity configuration.

Birefringent Filter—the micrometer adjustment protruding through the top cover tunes the wavelength by rotating the birefringent filter assembly.

Rod translation control—translates the Ti:sapphire rod in a plane parallel to the Brewster-angle face of the rod. By moving the rod in this manner, the cavity and pump modes avoid small flaws and bubbles within the rod. The rod can be translated to maximize the output power while the laser is operating. The translation stage must be released before attempting to move the rod. The stage is secured by two lock-down screws, one on each side (Figure 4-1).




Figure 4-1: Rod Translation Stage Lock-down Screws





Figure 4-2: Model 3900S Optical Element Identification



The horizontal control for mirror mounts M_3 to M_8 is the unlabeled knob.

Figure 4-3: Model 3900S Mirror Mount and Control Locations

Model 3910 Purge Unit Controls

The optional *Model 3910* contains oil and air filters and a molecular-sieve assembly to dry and clean nitrogen gas for purging the *Tsunami*[®] laser head. It has one control and two connectors. Refer to Figure 4-4.

FLOW ADJUST knob—used to set the laser head purge rate (gas flow) from 0.3 to $3.0 \text{ m}^3/\text{hr}$ (1.0 to 10.0 SCFH).

Flow adjust indicator—displays the nitrogen gas flow rate.

OUTPUT connector—connects the *Model 3910* to the laser head to provide dried, filtered, nitrogen gas.

INPUT FROM GAS CYLINDER connector—connects the *Model 3910* to a dry nitrogen gas supply. Limit the input pressure to 67 kPa (10 psi).




Figure 4-4: Model 3910 Regulator/Filter Purge Unit

Installation

Installing the Pump Laser

- Mount the pump laser on a suitable table.
 A honeycomb or granite optical table is best.
- 2. Set up the *Millennia*® series or argon ion pump laser according to its user's manual.
- 3. Adjust the feet of the pump laser so that its output beam is parallel with the table and is centered on the input widow of the *Model 3900S*. If there is not enough adjustment to lower your pump laser appropriately, raise the *Model 3900S*.
- 4. Verify the output of the pump laser meets specifications for power and mode quality.
- Reduce the pump laser output to a minimum.For ion lasers, reduce the tube current and close down the intracavity aperture.
- 6. If a large frame ion laser is used, use two mirrors to route the pump beam to the *Model 3900S* as shown in Figure 4-5. If a small *Millennia* laser is used, it can be set up so that it is in-line with the *Model 3900S*. In a typical ion laser setup, the routing mirrors will be about 15 cm in
 - In a typical ion laser setup, the routing mirrors will be about 15 cm in front of each laser and about 38 cm apart. We recommend using the two routing mirrors because it makes it much easier to direct the pump beam into the *Model 3900S*. Routing mirror holders and optics are available from Spectra-Physics. Refer to Service and Repair: Replacement Parts for specifications and part numbers.

Enclose the pump beam with beam tubes to avoid accidental exposure.




Figure 4-5: Typical system layout for *Model 3900S* and pump laser.

7. Close the pump laser shutter.

This completes the installation of the pump laser.

Installing the Model 3900S Ti:sapphire Laser

- 1. Place the Ti:sapphire laser on the table as shown in Figure 4-5.
- 2. Remove the top cover and raise the interlock shutter lever to the over-ride (up) position.
- 3. From inside the *Model 3900S*, use a ⁵/₃₂ in. Allen wrench to adjust the feet so the base plate of the laser is level with the table and the center of the input aperture is at the level of the pump beam.
- 4. Open the shutter of the pump laser.
- 5. Position the *Model 3900S* so that the pump beam passes unobstructed through the center of both the input aperture and the input alignment aperture, striking the approximate center of mirrors M₁, M₂ and M₃ (Figure 4-2).
- 6. Close the intracavity shutter on the pump laser.
- 7. Secure the *Model 3900S* to the table with foot clamps (provided). The clamps slide over the lower portion of each foot.
- 8. Attach two cooling lines to the fittings on the side of the laser (Figure 4-3), then attach one of the tubes to a filtered water supply and the remaining tube to a drain. Since the ion laser also requires a considerable water cooling supply, a small portion of this water may be diverted to the *Model 3900S* by placing "T" fittings in its supply and return lines.

Refer to Chapter 8, "Service and Repair: Replacement Parts," for the part numbers of appropriate fittings for your specific ion laser.

This completes the installation of the *Model 3900S* laser.

Installing the Optional Model 3910 Regulator/Filter Purge Unit

1. Verify the hoses inside the *Model 3910* are still properly connected to the plastic filters.

Turn the *Model 3910* upside-down and visually inspect for loose connections.

2. Place the *Model 3910* in a convenient place within 3.5 m of the *Model 3900S* laser head.

Warning!



PTFE tubing is included in the accessory kit. Use it when purging the *Tsunami* laser with dry nitrogen. Using PTFE tubing avoids introducing outgassed impurities into the cavity that may degrade system performance and/or cause damage to the optical coatings.

3. Connect the purge line from the *Model 3910* purge unit to the *Model 3900S* laser head.

Use the 3.5 m PTFE purge line provided to connect the *Model 3910* to the laser head. It comes with a Clip-Lok® connector on one end and a quick-disconnect connector on the other. The purge line is attached by pressing the hose connectors onto their mating connectors until they snap in place (Figure 4-4).

- a. Attach the Clip-Lok connector to the output of the *Model 3910*.
- b. Attach the quick-disconnect connector to the purge input on the *Model 3900S* laser head input bezel.

Warning!



Do not connect the *Model 3910* to a boil-off tank or other non-lab-grade nitrogen gas source, nor to any source with a pressure greater than 67 kPa (10 psi) or damage to the filters will result.

4. Connect the gas supply line from the *Model 3910* purge unit to the nitrogen supply. Make sure you use dry, oil-free, Electronic Grade 5 (or better) nitrogen (99.999% pure) to prevent contamination of the *Model 3910* or your laser optics.

Use the 2 m, ¹/₈ in. i.d. PTFE gas supply line provided to connect the *Model 3910* to the nitrogen supply. It comes with a Clip-Lok connector on both ends.

- a. Attach one end of the gas supply line to the INPUT FROM GAS CYL-INDER connector on the *Model 3910*.
- b. Attach the other end of the gas supply line to your nitrogen supply. If your nitrogen supply does not accept a Clip-Lok connector, cut off the existing connector and install your own. A Clip-Lok connector for ¼ in. i.d. PTFE tubing is also included in the accessory kit.

Clip-Lok® is a registered trademark of Anarak, Inc.

To attach barbed connectors to PTFE tubing, heat the tubing with a flameless heat source (heat gun), then quickly slide the tubing onto the barbed fitting while the tubing is still warm. Do not move the connection until it has cooled. Install the line and check for leaks.

This completes the installation of the *Model 3910* purge unit.

Removing the Purge Line from the Laser Head and Model 3910

To release the purge line from the laser head, press inward on the metal locking flange on the quick-disconnect, while gently pulling on the hose.

To release the purge line from the *Model 3910* purge unit, press inward on the wire clip-spring while gently pulling on the hose.

When you are ready to begin laser alignment, turn on the water supply and warm up the pump laser for 20 minutes.

Cavity Alignment

Note



Although there are many similarities between jet dye lasers and the Ti:sapphire laser, cavity alignment technique is not one of them. If you are familiar with dye laser alignment, do not attempt to apply these techniques to the *Model 3900S*. Please carefully follow the instructions below for best results.

Aligning the Mirrors and Rod

- 1. Place a power meter in the beam path at the output of the Ti:sapphire laser. Verify that the pump laser is set to minimum power.
- 2. Locate initial optics positions.

Note



Although all optics positions were optimized at the factory, all alignment procedures are repeated here in the event that misalignment has occurred.

Release the lock-down screws on the rod translation stage. Move the rod back and forth using the translation control, noting the extent of allowed travel. Adjust the rod to the center of this range and gently secure.

Warning!



The alignment of both the birefringent filter and Ti:sapphire rod have been carefully optimized at the factory. Under no circumstances should either be removed from their mounts, or their alignments altered, except for a translation adjustment of the rod. Any other modification or adjustment made by persons other than a qualified Spectra-Physics service engineer will void the warranty.

Release the lock-down screws securing fold mirrors M_5 and M_6 so they can be adjusted by gently sliding the mount (Figure 4-6). The initial position of each mount should be as shown; the edges of the fixed and adjustable portions of the slide mounts are offset by about 0.5 mm. The correct distance from the centers of M_5 and M_6 to their respective nearest rod face is 47 mm. Lock the mirror slides in position once the adjustment is complete.




Figure 4-6: Location of M₅ and M₆ lock-down screws.




Figure 4-7: M₄ Lock-down Screws and Micrometer Setting

Release the lock-down screws securing the pump focus mirror M_4 . Adjust the micrometer on the M_4 mirror mount so the vernier reads 3.5 (Figure 4-7). This will set the distance from the front face of M_4 to the rear face of M_5 to approximately 69 mm. The optimum pump focus mirror position will depend somewhat on the divergence of the pump beam and the distance between the pump and Ti:sapphire lasers. Do not secure the pump focus mirror at this time.

3. Open the pump laser shutter.

Verify that the pump beam still strikes the center of M_1 and M_2 and that the pump power is at a minimum.

4. Align M₃ and M₄.

Orient M₃ horizontally (green) and vertically (blue) so the pump beam strikes the approximate center of the face of M₄

Orient M₄ horizontally and vertically so the focused pump beam passes through the centers of the concave surfaces of M₅ and M₆, striking the center of the Ti:sapphire laser rod. The concave surfaces are those facing the Ti:sapphire rod. You might need to adjust the translation of the Ti:sapphire rod.

Note



By this point, the pump beam will have interacted with a total of nine coated and polished surfaces, including the Ti:sapphire rod itself. Dirt or dust particles on these surfaces will degrade the pump beam. If you suspect dirt on any surface, refer to Chapter 7, "Maintenance," for optics cleaning instructions.

5. Raise the power of the pump laser.

For safety, close both Ti:sapphire and pump laser shutters. Raise the pump laser tube current to approximately 4 or 5 W and let it warm up for 10 minutes before making the remaining adjustments. During this time, ensure that the faces of the Ti:sapphire rod are clean. Finally, open both shutters.

6. Adjust M₅.

Place a white card over the aperture at M₇ and adjust M₅ to center a broad fluorescence spot on M₇ as shown in Figure 4-8.

7. Adjust M_6 .

Place the white card over the aperture of M_8 Some of the pump beam will be reflected from M_6 toward M_8 Although the pump beam will not follow exactly the same path as the cavity mode, this small blue spot can be used as a reference. Center this spot on the aperture of M_8 by adjusting M_6 .




Figure 4-8: The fluorescence spot as it appears on mirror M₇.

Note



A pair of argon goggles (laser goggles or glasses opaque to argon laser wavelengths) is an ideal tool for this adjustment as this will make the residual pump beam appear faint while making the fluorescence from the Ti:sapphire more obvious. The fluorescence from the Ti:sapphire rod will appear more diffuse at M_8 than it will at M_7 , but the adjustment technique is the same as it was for M_5 .

8. Set the birefringent filter.

Adjust the operating wavelength of the birefringent filter to correspond to 800 nm. If the orientation of the filter has not been changed since factory installation, this position should correspond to about "4.8" on its micrometer vernier (Figure 4-3). The blue dot on the outside edge of the filter should align with the setscrew on the filter adjustment.

9. Adjust M₈.

Place the white card over the aperture to M_7 . Adjust M_8 to center a small fluorescent dot on the aperture. This fluorescent spot can be seen with the naked eye by many people. However, if the pump power is rather low or your visual acuity is poor in the extreme red, try using the argon goggles described in Step 7, or use an infrared viewer.

10. Adjust M₇.

Place the white card over M₈ and adjust M₇ to center a small fluorescent spot (fainter than that in Step 9) onto M₈.

Remove the card and increase pump power to the recommended operating power; the *Model 3900S* should lase. If not, place the white card over the power meter and observe the faint spots on the card with the goggles or infrared viewer. A pattern of misaligned spots should be visible that can be aligned by adjusting M₈ and M₇ to start the system lasing.

11. Adjust for maximum power output: M₇ and M₈.

Before this centering is attempted, adjust the laser for maximum power. To achieve maximum power output, M_7 and M_8 must be "walked" and the laser mode must strike the center of M_7 and M_8 .

Adjust only M_7 and M_8 to maximize power output. They must be adjusted in an organized manner. Adjust the horizontal (green) and vertical (blue) controls of M_7 and M_8 in pairs.

First, adjust the horizontal control of M_8 for maximum power. Then adjust the horizontal control for M_7 for maximum power. Repeat this procedure back and forth until no further increase in power is possible from horizontal adjustments alone.

In a similar manner, adjust only the vertical controls of M_7 and M_8 for maximum power.

Go back and forth between the horizontal and vertical controls in pairs until maximum power is achieved.

12. Adjust for beam centering.

Verify that the intracavity beam is centered on M₇ and M₈ and that it is unobstructed. If so, skip to Step 13.

To center the beam on M_7 or M_8 , or both, the fold mirrors will have to be adjusted. However, adjust only one fold mirror at a time: M_5 and M_7 will be adjusted as a pair, as will M_6 and M_8 . This technique leaves one of the two collimated legs of the four-mirror cavity unaffected, making alignment much easier.

To center the beam on M_7 , slightly adjust M_5 to move the intracavity beam in the desired direction. The unit will stop lasing, but can be re stored by adjusting M_7 After the beam is centered and lasing is re stored, maximize power by walking M_7 and M_8

Repeat this last step for M₆ and M₈

13. Adjust the rod position.

To maximize output power in the rod, carefully release the rod translation lock-down screws and slowly translate the rod back and forth while monitoring the output power. Select a rod position that yields maximum output, then gently tighten the translation stage lock-down screws.

Scattering centers in the rod can often be located by observing the residual pump beam that is normally deposited in the beam dump. Temporarily close the interlock shutter and remove the entire beam dump from the baseplate. Place a screen to intercept this residual pump beam several meters from the laser and open the shutter. When the rod is translated, flaws and bubbles will appear as non uniformities in the spot.

14. Adjusting for maximum power: M₄ position.

If M_4 has not been positioned according to factory specification, the laser may display **two** points where maximum power is obtained when horizontally scanning M_7 and M_8 (Step 11 above).

Note



The pump mirror (M_4) should only be moved within limits that allow the laser to continue to lase. This may be several turns of the micrometer adjustment if the laser is well above threshold, but only one turn should be attempted at a time.

To adjust, note the laser output power, then turn the micrometer to a new position and carefully adjust mirrors M_7 and M_8 as discussed in Step 11. If more power is obtained, continue with further adjustment of the pump focus in the same direction. If power is reduced, reverse the direction of adjustment. Remember: adjust by no more than one full turn at a time.

Adjusting the Birefringent Filter

Warning!



The birefringent filter assembly is carefully aligned at the factory to orient the filter plates at Brewster's angle. Never remove this assembly from the baseplate. Also, never attempt to disassemble the filter stack.

The birefringent filter has both a blue dot and a yellow dot on its side. For operation between 700 to 850 nm, align the blue dot with the setscrew in the side of the filter adjustment. Reference the yellow dot for operation between 850 to 1000 nm. Loosen the setscrew to rotate the birefringent filter within the mount from one dot to the other. Take care not to over tighten the setscrew.

Observe the beam through a reversion spectroscope or a monochromator and infrared viewer. Set the micrometer adjustment to "7.5" when using the 700 to 850 nm mirror set or to "0.5" when using the 850 to 1000 nm set (refer to Figure 4-9). With the spectroscope or monochromator set to 850 nm, adjust the orientation of the filter until the laser operates at 850 nm, then tighten the setscrew. The laser will now tune approximately according to the calibration chart in Figure 4-9. For accurate measurement of wavelength, use a monochromator to fine tune the laser to the desired operating wavelength.

Optional sets of optics can be purchased to broaden the tuning range of the *Model 3900S* system. Table 8-1 and Table 8-2 give a complete listing of all the optic sets available for this laser along with the birefringent filter setting. Align the appropriate colored dot on the side of the birefringent filter as outlined above under "Adjusting the Birefringent Filter."

A different birefringent filter is required for use with the 800 to 900 nm optic set, as well as for use when near or beyond 1100 nm (with large frame pumping). This filter has purple and white dots on the side of the assembly, giving a reference rotation point.

The calibration curve (Figure 4-10) gives an approximate wavelength vs. micrometer setting for the various optional optic sets. However, a wave length measuring device should be used to obtain an accurate reading.



Figure 4-9: Wavelength Calibration Chart




Figure 4-10: Optional Wavelength Calibration Chart

Two different solid etalons may be added to the *Model 3900S* Ti:sapphire cw tunable laser to achieve narrow linewidth operation. The use of two equal length collimated legs with a centered gain medium allows laser operation in two adjacent longitudinal modes with little power loss.

The thin etalon (0.5 mm thick with a free spectral range (FSR) of 200 GHz) allows operation with less than 15 GHz linewidth. It is tuned by tilting the etalon via a black knob on the rear panel. Continuous tuning is possible over a range larger than the nominal bandwidth of the birefringent filter (approximately 100 GHz).

The thick etalon (5 mm thick with a FSR of 20 GHz) allows operation in two or three adjacent cavity modes spaced by about 200 MHz. Tuning is accomplished by either angular adjustment of the etalon or by temperature control of the etalon housing.

Thin Etalon Assembly

Installation

Figure 5-1 and Figure 5-2 show the location of both the thick and thin etalon assemblies. Figure 5-3 shows how mirror mount M₈ is modified to accept the thin etalon mount.

1. Optimize laser performance.

Before installing any etalons, optimize the laser performance as described in Chapter 4, "Installation and Alignment: Cavity Alignment."

2. Remove the adjustment knobs for M₈.

Close the shutter on the pump laser. Pry off the caps on both knobs for the mirror mount M₈. Remove the knobs by loosening their collets with a screwdriver.

3. Remove the rear escutcheon.

Remove the rear escutcheon along with the attached shutter assembly and input beam trim ring (four screws). Retain these screws. If beam tubes were used in the installation, slide the final beam tube into the trim ring to allow the escutcheon to be removed without disturbing any turning mirrors.

4. Remove the M₈ mirror mount.

Remove the entire M₈ mirror mount by loosening screws marked "A" in Figure 5-3. These screws will be used to install the etalon mount.

5. Remove the adapter plate.

Loosen screws marked "B" and remove the adapter plate from the M₈ mount.

6. Assemble the etalon mount and M₈ mount.

Use screws marked "B" to attach the M₈ mirror mount to the etalon mount assembly. Use screws marked "A" to install the assembly on the baseplate.

7. Assemble the escutcheon.

Remove the shutter assembly and trim ring from the old escutcheon and install them on the new plate, then install on the rear bezel.

8. Reinstall the knobs and knob caps for M₈.

Alignment

1. Realign laser.

The etalon should not yet be installed into its mount. Since M_8 was removed to assemble the etalon mount, realign it according to Step 9 of "Cavity Alignment: Aligning the Mirrors and Rod," in Chapter 4. Adjust M_8 for maximum output power.

2. Set up a scanning interferometer.

To fully utilize the thin etalon, measure the linewidth of the laser output using a scanning interferometer, such as a Spectra-Physics *Model 410* scanning etalon or equivalent, with a FSR of about 300 GHz and a linewidth of approximately 5 GHz. Use a thick, wedged beam splitter to reflect a portion of the output into the interferometer and optimize according to its instruction manual.

Clean etalon.


Make sure the etalon is cleaned thoroughly before installing. Use lens tissue and methanol in the same manner as used to clean the laser mirrors.

4. Install the etalon.

Install the etalon into its mount. If the etalon is perpendicular to the intracavity beam, the unit will lase when the pump laser shutter is opened.

5. Adjust the etalon tilt.

Slightly release the knurled nut near the vertical adjustment knob and, using an Allen wrench, adjust the etalon vertically for maximum output power. Tighten the knurled nut to secure.



The horizontal control for mirror mounts M_3 to M_8 is the unlabeled knob.

Figure 5-1: *Model 3900S* Mirror Mount and Control Locations, including Etalon Assembly Options




Figure 5-2: Model 3900S with Etalon Assembly Options




Figure 5-3: Etalon Mount Assembly with M₈ Modification

6. Tune the laser output.

Adjust the black knob on the rear escutcheon to tune the laser.

Bear in mind the following as you use the thin etalon:

- Cavity alignment can change the cavity mode beam direction with respect to the thin etalon. This, then, will change the wavelength of the laser, even though the birefringent filter has not been adjusted.
- Although the linewidth obtained with the birefringent filter is less than 40 GHz, the etalon will typically tune over a range of 100 GHz.
- An etalon placed perpendicular to the beam will cause instability in the operating wavelength as evidenced by jitter in the scanning interferometer trace. Slightly readjusting the vertical tilt adjustment knob as per Step 5 above will eliminate this instability.

This completes the installation and alignment of the thin etalon assembly.

Thick Etalon Assembly

Installation

The thick etalon is used with the thin etalon to achieve sub-gigahertz linewidth from the *Model 3900S*. Figure 5-1 and Figure 5-2 show the location of both the thick and thin etalon assemblies.

1. Optimize laser performance.

Before installing the thick etalon, optimize the laser performance.

2. Set up a scanning interferometer.

To fully utilize the thick etalon, use a scanning interferometer, such as a Spectra-Physics *Model 470* or equivalent, with a FSR of about 2 GHz and a linewidth of approximately 20 MHz. Use a wedged beam splitter to reflect a portion of the output into the interferometer; optimize the alignment according to the interferometer manual.

3. Remove the adjustment knobs for M₈.

Close the shutter on the pump laser. Pry off the caps from both knobs with a thin blade. Remove the knobs by loosening their collets using a screwdriver. If done with care, this mirror will need minimal realignment.

4. Remove the rear escutcheon.

Remove the rear escutcheon, the shutter assembly, and input beam trim ring (4 screws). If beam tubes have been used in the laser installation, slide the final beam tube into the trim ring to allow the escutcheon to be removed without disturbing any turning mirrors.

5. Assemble the new escutcheon.

Remove the shutter assembly and trim ring from the old escutcheon and install on the plate provided with the etalon. Install onto the rear bezel.

- 6. Install the knobs and knob caps for M₈.
- 7. Align the laser.

Before installing the thick etalon, align M_8 only according to Step 9 of "Cavity Alignment: Aligning the Mirrors and Rod," in Chapter 4. Open the pump laser shutter and align M_8 for maximum output power.

8. Install the thick etalon.

Close the pump laser shutter and install the thick etalon assembly to the base plate (2 screws, Figure 5-1).

9. Connect electrical cables.

Connect the etalon cable to the one from the rear escutcheon. Connect the receptacle on the escutcheon to the *Model 482* etalon temperature controller. Plug the power cord into the *Model 482* and to your power source.

Alignment

1. Turn on the temperature controller.

Adjust the temperature control to "5." Let the etalon warm up until the panel light turns off.

2. Optimize power output.

If the thick etalon is reasonably perpendicular to the intracavity beam, the unit will lase when the pump laser is turned on. However, because of the finesse and thickness of the etalon, its orientation may radically affect output power. Adjust the etalon vertically and horizontally for maximum output. Then turn the blue knob about a ¼ turn away from the maximum power position. This prevents wavelength instability caused by feedback into the laser cavity.

3. Optimize mode structure.

Use the scanning interferometer to observe the linewidth. When operating perfectly, the interferometer will indicate pairs of longitudinal modes separated by the FSR of the interferometer (2 GHz) with a mode spacing of 200 MHz. If the laser cavity is misaligned, multiple transverse modes will be exhibited as small "bumps" alongside the main longitudinal modes. Further, as pump power is increased, another longitudinal mode will often appear. Tilt the etalon to reduce this second mode.

If the etalon is exactly perpendicular to the intracavity beam, the mode structure will be unstable. To eliminate this instability, slightly misalign the etalon vertically, horizontally, or both.

4. Tune the laser output.

Because tilting the etalon affects mode stability, it is best to temperature tune the output. Increasing the temperature of the etalon length ens the solid etalon, decreasing the operating frequency. After adjusting the temperature, wait until the temperature has stabilized before expecting the operating wavelength to remain fixed. Wavelength change may be monitored with the scanning interferometer.

This completes the installation and alignment of the thick etalon assembly.

Chapter 6 Operation

System Startup

Once the *Model 3900S* has been installed and aligned, the day-to-day operation is consistent and repeatable. Provided nothing is disturbed between operating sessions, only minor adjustments are required to return the system to optimum performance levels.

- 1. Turn on the cooling water supply.
- 2. Close the shutter on the pump laser, then turn on the laser at low cur rent and allow it to warm up. Refer to the ion laser manual for the time required for warm up.
- 3. Increase the current to the anticipated operating level and allow the laser to warm up for an additional 5 minutes.
- 4. Begin the dry nitrogen purge.

The optional *Model 3910* regulator/filter purge unit is for use with bottled nitrogen gas to eliminate problems associated with dust, contamination and the tuning discontinuities caused by oxygen and water vapor. The latter exhibits strong absorption lines in the long wave length region between 840 and 1000 nm, especially at 944, 950, 981.2 and 984.1 nm (Figure 6-1). Due to the strength of these absorption lines, the moisture level within the laser head must be below 1000 ppm to allow problem-free mode locking. For these wavelengths, we recommend using 99.999% pure, dry, oil-free, Electronic Grade 5 nitrogen at a starting purge rate of 0.3 m³/hr (10 SCFH) for at least 2 hours, then reducing it to 0.17 m³/hr (6 SCFH) for the duration of the experiment. The purge can be reduced to 0.014 m³/hr (0.5 SCFH) for all other wavelengths.

- a. Verify the laser head cover is in place and clamped down, sealing the cavity.
- b. Verify the nitrogen tank output regulator is set to minimum, then turn on the nitrogen supply.
- c. Set the output regulator to limit pressure to less than 67 kPa (10 psi).
- d. Open the *Model 3900S* purge bleed valve.
- e. Use the *Model 3910* flow control to set the nitrogen flow rate as outlined above.




Figure 6-1: Transmittance vs. Wavelength for Oxygen and Water Vapor.

- 5. Open the pump laser shutter and note the output. Optimize the pump laser power if necessary.
- 6. With a power meter in the output beam of the *Model 3900S*, adjust the following mirrors for optimum performance:
 - a. Adjust M₈ and M₇ horizontally until power is maximized.
 - b. Adjust M₈ and M₇ vertically until power is again maximized.
 - c. Repeat these alternating horizontal and vertical adjustments until no further power increase is observed.

System Shutdown

- 1. Turn off the pump laser, following the shutdown procedure in its instruction manual.
- 2. Turn off the cooling water. To prevent condensation of the surface of the rod, do not circulate cooling water when the laser is off.

Tuning with the Birefringent Filter

Use an infrared viewer and a monochromator when tuning the *Model 3900S*. Center the output of the *Model 3900S* on the input slit of the monochromator and set the wavelength of the monochromator to the desired operating wavelength. Use the micrometer adjustment to tune the laser to within the resolution of the monochromator. If you wish to change the operating wavelength range, you will need to change both the optics set and the order of the birefringent filter.

Optics Selection

Five mirror sets and wavelength operating ranges are available for the *Model 3900S*. A complete list of the mirror options available for the *Model 3900S* can be found in Table 8-1 and Table 8-2 in Chapter 8, "Service and Repair: Replacement Parts."

Changing Mirrors

Two sets of optics allow the laser to operate over wavelengths from 700 to 850 nm or 850 to 1000 nm. It is easiest to switch optics after the laser is operating and optimized. When changing the wavelength range, change the fold mirrors (M_5 and M_6), the output mirror (M_7), and the high reflector (M_8) one at a time, in order. Close the pump laser shutter before replacing each optic and open the shutter only after each is cleaned and installed.

- 1. Tune the laser to 850 nm and optimize the output power.
- 2. Replace the mirrors one at a time and optimize the laser after instal ling each. Clean each mirror before installing. First replace M₅. Realign M₇ to maximize power.
- 3. Next replace M₆ Only adjust M₈ to restore power.
- 4. Next replace M₇ and adjust for maximum power.
- 5. Finally, replace M₈ and adjust for maximum output power. At this point, you may wish to optimize output power as outlined in Chapter 4, "Installation and Alignment: Cavity Alignment."
- 6. When changing optics, the filter order must be changed. Rotate the filter in the assembly, aligning the appropriate colored dot on the filter with the setscrew on the filter adjustment assembly: the blue dot for 700 to 850 nm, the yellow dot for 850 to 1000 nm. Turn the micro meter to "7.5" for 700 to 850 nm or to "0.5" for 850 to 1000 nm. Using a monochromator and an infrared viewer, rotate the filter in the adjustment assembly until the laser operates at 850 nm. At this point, the tuning curve will closely follow the wavelength calibration shown in Figure 3-9. Tighten the setscrew until barely snug. Adjust M₇ to optimize output power after changing the filter order.

Instructions for removing and cleaning mirrors are found in Chapter 7, "Maintenance."

Chapter 7 Maintenance

The condition of the environment and the amount of time the laser is used will affect your periodic maintenance schedule. Optics will obviously stay clean much longer if not exposed to smoke or other air-born contaminants. Condensation due to excessive humidity can also contaminate the optical surfaces. Try to provide a smoke-free, filtered, dry environment for the laser. The cleaner the environment, the slower the rate of contamination.

Notes on the Cleaning of Laser Optics

Ion lasers are oscillators that operate with gain margins of a few percent. Losses due to unclean optics, which might be negligible in ordinary optical systems, can disable a laser. Dust on mirror surfaces can reduce output power or cause total failure. Cleanliness is, therefore, essential. The maintenance techniques used with laser optics must be applied with extreme care and with attention to detail.

"Clean" is a relative description; nothing is ever perfectly clean, and no cleaning operation ever completely removes contaminants. Cleaning is a process of reducing objectionable materials to acceptable levels.

Since cleaning simply dilutes contamination to the limit set by solvent impurities, solvents must be as pure as possible and leave as little solvent on the surface as possible. As any solvent evaporates, it leaves impurities be hind in proportion to its volume. Avoid rewiping a surface with the same lens tissue; a used tissue and solvent will redistribute contamination, they won't remove it.

Always use fresh solvent. Both methanol and acetone collect moisture during prolonged exposure to air. Avoid storage in bottles where a large volume of air is trapped above the solvent. Instead, store solvents in small glass bottles where either the solvents are used up quickly or the bottles are filled frequently from a fresh, uncontaminated source.

Laser optics are made by vacuum depositing microthin layers of materials of varying indices of refraction on glass substrates. If the surface is scratched to a depth as shallow as 0.01 mm, the operating efficiency of the optical coating will be reduced significantly.

Because an intracavity passive catalyst is used in the Chroma 5 laser, there should be little requirement for optical cleaning if the cavity is left closed and undisturbed. However, because optics are exchanged from time to time, there will be occasions when cleaning is necessary.

Stick to the following principles whenever you clean any optical surface.

- Remove and clean one optic at a time, then replace it and optimize laser output power. If more than one optic is removed at time, all reference points will be lost, making realignment extremely difficult.
- Perform any maintenance in a clean environment, over an area covered by a soft cloth or pad, if possible.
- Wash your hands thoroughly with liquid detergent. Body oils and contaminants can render otherwise fastidious cleaning practices useless.
- Use dry nitrogen, canned air, or a rubber squeeze bulb to blow dust and lint from the optic surface before cleaning it with solvent. Permanent damage may occur if dust scratches the coating.
- Use spectrophotometric-grade (HPLC) solvents. Do not try to remove contamination with a cleaning solvent that may leave other impurities behind.
- Use powder-free, clean latex gloves or finger cots.
- Use Kodak Lens Cleaning Paper[™] or equivalent to clean optics and plasma tube windows. Use each piece only once: a dirty tissue merely redistributes contamination, it does not remove it.

Warning!



Do not use lens tissue that is designated for cleaning eye glasses. Such tissue contains silicones. These molecules bind themselves to the mirror coatings and window quartz and can cause permanent damage. Also, do not use cotton swabs, e.g., Q-TipsTM. Solvents dissolve the glue that is used to fasten the cotton to the stick, and the result is contaminated coatings. Only use photographic lens tissue to clean optical components.

Equipment Required

Some of the following are supplied in the accessory kit.

- Dry nitrogen, canned air, or rubber squeeze bulb
- Plastic hemostat
- Clean (new) finger cots or powder-free latex gloves
- Kodak Lens Cleaning Paper or equivalent





Do not use a metal hemostat or forceps for cleaning optics or windows. There is danger of electrocution when using these instruments in side the laser head. Use the plastic hemostat provided.

[&]quot;Kodak Lens Cleaning Paper" is a product of the Kodak Corporation "Q-Tips" is a product of the Johnson and Johnson Corporation

Cleaning Solutions Required

• spectrophotometric-grade (HPLC) acetone or methanol

Warning!



Only use spectrophotometric-grade acetone and/or methanol. Using lower grade reagents may lead to contamination of or damage to optical coatings.

Note



Methanol tends to clean better but, if not fresh, may deposit a water-based film on the surface being cleaned. If this occurs, follow the methanol wipe with an acetone wipe to remove the film. As always, use fresh solvent from a bottle with little air in it.

General Procedures for Cleaning Optics

The laser should be on and stable. Remove, clean, and replace one optic at a time, maximizing laser output power after each optic is cleaned. Use clean finger cots to protect all intracavity components, including the coated mirror surface.

Do not allow the cavity to remain open for very long. The intracavity passive catalyst that reduces O_3 in the cavity will become contaminated and, therefore, ineffective. Keep the output coupler and high reflector locked in place, and make sure the tubular cavity seals are in place. Following these simple rules will ensure a long catalyst life.

The high reflector can be left in its holder for cleaning. However, the output coupler (OC) and the beam splitter must be removed from their holders to clean the second surface.

1. Use a squeeze bulb, dry nitrogen, or canned air to clean away any dust or grit before cleaning optics with solvent.

Drop and Drag

- 2. Whenever possible, clean the optic using the "drop and drag" method (Figure 7-1).
 - a. Hold the optic horizontal with its coated surface up. Place a sheet of lens tissue over it and squeeze a drop or two of acetone or methanol onto it.
 - b. Slowly draw the tissue across the surface to remove dissolved contaminants and to dry the surface.
 - Pull the tissue slow enough so the solvent evaporation front immediately follows the tissue, i.e., the solvent dries only after leaving the optic surface.



Figure 7-1: Cleaning the Mirror Surface

Tissue in Hemostat

- 3. For stubborn contaminants and to access hard-to-reach places (such as the windows), use a tissue in a hemostat to clean the optic.
 - a. Fold a piece of tissue in half repeatedly until you have a pad about 1 cm (0.5 in.) square, and clamp it in a plastic hemostat (Figure 7-2).

Don't Touch!



While folding, *do not touch* the surface of the tissue that will contact the optic, or you will contaminate the solvent.




Figure 7-2: Lens Tissue Folded for Cleaning

- b. If required, cut the paper with a solvent-cleaned tool to allow access to the optic.
- c. Saturate the tissue with acetone or methanol, shake off the excess, resaturate, and shake again.
- d. Wipe the surface in a single motion.
 Be careful that the hemostat does not touch the optic surface, or the coating may be scratched.
- 4. After placing the optic you just cleaned into the beam, inspect it to verify the optic actually got cleaner, i.e., you did not replace one contaminant with another.

Steps for Cleaning Optics

- 1. If necessary, remove the optic to be cleaned from its holder to gain access to optic surface(s).
- 2. For easy access optics, hold the optic horizontal with its coated surface up, place a sheet of lens tissue over it, and squeeze a drop or two of acetone or methanol onto it. Gently draw the tissue across the surface to remove dissolved contaminants and to dry the surface.
- 3. If stubborn contaminants remain, use the "Tissue in Hemostat" procedure described above to clean the surface.
- 4. If the optic is an output coupler, invert the optic and repeat the procedure to clean the second surface.
- 5. For the Brewster windows, use the "Tissue in Hemostat" procedure, but cut the tissue to fit within the Brewster window slot using clean scissors or wire cutters.

Removing the Mirrors for Cleaning

Remove, clean and install mirrors one at a time to avoid accidental exchanges and misalignment. Close the pump laser shutter while optics are being changed. Further, adjust the laser for maximum output after cleaning each mirror.

Output Mirror M₇

Close the pump laser shutter, then unscrew the knurled mirror holder to remove the optic holder. You need not remove the mirror from its holder to clean the coated surface.

If the front surface of the mirror needs cleaning, grasp the edges of the mirror and gently work it out of the holder; use finger cots to prevent contamination. After cleaning, screw in the mirror holder until snug; finger tight is sufficient. Excessive tightening will damage the optic.

After the optic is installed, open the pump laser shutter and adjust only M₇ to restore maximum power.

High Reflector M₈

Clean and replace the optic in the high reflector mirror holder in the same fashion as described for M₇. After installation, adjust only M₈ for maximum power.

Fold Mirrors M₅ and M₆

A knurled, spring-loaded sleeve holds each of these mirrors in place. Unscrew the sleeve from the mirror holder and clean each optic within its holder. After cleaning, carefully install the mirror, making sure that the mirror is not binding and is actually free to move side to side slightly, as it rests against the three balls within the mirror mount. As with M₇ and M₈, screw the optic in place until snug. These mirrors must be cleaned with

great care to insure optimum laser output. After installing M_5 , open the pump shutter and adjust M_7 for maximum power. Better control over alignment is achieved through an adjustment of M_7 rather than M_5 , unless the orientation of the mirror mount for M_5 has been disturbed.

In the same fashion, clean and replace M₆ next, and adjust M₈ for maxi mum output.

Pump Mirrors M_1 , M_2 , M_3 and M_4

This mirror set is not changed when wavelength ranges are changed and, since they are not part of the laser cavity and are less sensitive to contamination, they will not often need to be removed. Furthermore, since the pump beam alignment is very dependent on the orientation of these mirrors, removal and replacement of these mirrors may require more extensive realignment of the entire laser. Nevertheless, if these mirrors become excessively dusty, cleaning them may improve laser performance.

Cleaning the Ti:sapphire Laser Rod Surfaces

Laser performance is extremely sensitive to dirt on the rod surfaces. Do not remove the rod from the mount; a small portion of rod protrudes from the end of the holder, allowing easy access for cleaning.

After removing the beam shield from the top of the rod holder, clean the rod surface with a pad of lens tissue held in a hemostat and lightly saturated with methanol or acetone. Make only one pass over the rod surface each tissue; wipe in one direction only. Install the beam shield over the Ti:sapphire rod.

Cleaning the Birefringent Filter

Warning!



The birefringent filter assembly is carefully aligned at the factory to orient the filter plates at Brewster's angle. Never remove this assembly from the baseplate. The filter stack should never be disassembled to at tempt to clean the inner surfaces.

Typically, only the upward-facing surface of the birefringent filter tends to get dusty, and it may be cleaned without removing the filter from the system.

If the lower surface requires cleaning, remove only the filter stack from its mount. Note which color dot is aligned with the Allen screw in the mount, loosen the screw using the wrench provided in the laser tool kit, and gently pull the filter out.

Clean the surface(s) with a folded, acetone-saturated tissue clamped in a hemostat. Use very little pressure; the filter is made of thin material and can be very easily damaged.

If you removed the filter from its holder, make sure to install it with the appropriate dot aligned with the setscrew. Do not overtighten screw. If the

filter micrometer adjustment has been calibrated for wavelength, recalibration will now be necessary. Additionally, if the filter is removed, adjust M₇ for maximum power after installing the filter.

Cleaning the Brewster Windows

The *Model 3900S* incorporates Brewster angle windows at the input and output ends to seal the unit, allowing the system to be nitrogen purged. (Refer to Table 8-3 in Chapter 8 for purge options.)

Close the pump laser shutter. Remove the retaining screws on both sides of the trim rings, front and rear, and remove the inner sleeve to expose the Brewster windows for cleaning.

Clean the surfaces with a folded, acetone-saturated tissue clamped in a hemostat. When cleaning the inside surfaces, cut the excess folded tissue to allow the tissue and hemostat to fit through the window access slot in the holder.

Replacing Filters in the Optional Model 3910 Purge Unit

The schedule for replacing the filters and dryer/sieve assembly depends on the amount and quality of purge nitrogen used. Change all three filters (Figure 7-3) when the blue desiccant in the sieve assembly turns pink. In some areas, the indicator dye in the desiccant is considered a hazardous material. Consult your materials manager and/or local government environmental agency for guidelines on proper disposal methods. Refer to the Material Safety Data Sheets in Appendix A for information on the materials contained in these filters.

Part numbers for filters are listed in Table 8-3 at the end of Chapter 8.

To replace the filters and sieve assembly:

- 1. Turn the unit upside down with the flow gauge facing away from you, and remove the small screws (2) on the side nearest you.
- 2. Disconnect the filter assembly from the flow gauge and output port. Press in on the connector spring-clips and pull out the hoses.
- 3. Rotate the entire filter assembly so the small filters point upward, then move the assembly to the right (away from the flow gauge) and lift it out.
- 4. Remove the filter hose fittings.
 - Note the placement and orientation of each filter before disassembling the hoses. Loosen the screws and remove the filter assemblies.
- 5. Discard the filter assemblies.
- 6. Assemble a new filter assembly from new parts.

 The screws at each end of the two small filters should be tight. However, over-tightening will crack the plastic.
- 7. Place the new filter assembly in the box and connect and secure it in reverse order of disassembly.

Push the hose fittings into their mating connectors until they snap into place.

This completes the filter replacement procedure.




Figure 7-3: Underside of the *Model 3910*, showing filter assembly placement and orientation.

Troubleshooting Guide

Use this guide if *Model 3900S* performance drops unexpectedly or if performance cannot be restored after moving the laser or changing optics. If, after trying the following troubleshooting suggestions, you are still unsuccessful in returning your laser to normal operation, call your Spectra-Physics representative for help.

Symptom: No Output Power	
Possible Causes	Corrective Action
Insufficient pump power	Pump power must be at least 5 W.
Birefringent filter order is set incorrectly	Verify the birefringent filter is set to the correct order for the operating wavelength. Align the blue dot with the setscrew when using 700 – 850 nm optics.
Cavity optics are out of alignment	Adjust M_{7} and M_{8} as described in Chapter 4, "Installation and Alignment: Cavity Alignment."
Pump beam optics are out of alignment	This requires performing the complete alignment procedure as described in Chapter 4, "Installation and Alignment: Cavity Alignment."

Symptom: Low Output Po	wer
Possible Causes	Corrective Action
Low pump power	Check the pump laser output power and TEM ₀₀ mode.
Optics are dirty	Clean optics using the procedure described in Maintenance in the following order: 1. Rod surfaces 2. Birefringent filter, outside surfaces only 3. M ₈ 4. M ₇ 5. M ₆ 6. M ₅ 7. Brewster windows 8. M ₁ and M ₂ without removing the optic 9. As a last resort, remove and clean M ₃ and M ₄ . This will necessitate a full realignment. Adjust M ₇ and M ₈ for maximum output power after each stage of cleaning.

Symptom: Low Output Power		
Possible Causes	Corrective Action	
Cavity optics are out of alignment	Adjust mirrors M_7 and M_8 as described in Chapter 4, "Installation and Alignment: Cavity Alignment."	
Rod position are not optimized	Adjust the rod position as described in Chapter 4, "Installation and Alignment: Cavity Alignment." Do not forget to release the lock-down screws before adjustment.	
Pump focus mirror M ₄ out of position	Unless the position of M_4 has been moved since installation, this adjustment should not be necessary. However, if the Model 3900S is used with a pump laser that has a beam divergence differing greatly from the recommended Spectra-Physics argon ion lasers, M_4 can be adjusted as described in Chapter 4, "Installation and Alignment: Cavity Alignment."	

Replacement Parts

Table 8-1: Replacement Parts, Standard Model 3900S Optics

Optic	Description	Coating Range	Part Number
M_1, M_2, M_3	Pump beam flat	450-530 nm	G0342-001
M_4	Pump beam focus	475-530 nm	G0323-001
M_5, M_6	Cavity fold	700-850 nm	G0079-012
M_5, M_6	Cavity fold	850-1000 nm	G0079-013
M_7	Output coupler	700-850 nm	G0324-003
M_7	Output coupler	850-1000 nm	G0324-005
M_8	High Reflector	700-850 nm	G0324-002
M_8	High Reflector	850-1000 nm	G0324-004
Birefringent filter	Three-plate		0440-5530
Brewster window assembly			0441-8110

Table 8-2: Model 3900S Optic Sets

Wavelength	Mirror Set	Description	Birefringent Filter	
Range	Part No.		Part No.	Comments
675-750 nm ¹	0440-5505	Optional short mirror set	0434-8973	Align on yellow dot
700-850 nm	0440-5501	Standard short mirror set	0434-8973	Align on blue dot
800-900 nm ¹	0440-5503	Optional middle mirror set	0434-8974	Align on white dot
850-1000 nm	0440-5502	Standard long mirror set	0434-8973	Align on yellow dot
950-1100 nm ¹	0440-5504	Optional long mirror set	0434-8973 ²	Align on blue dot

Specifications only available after factory test (at additional charge) on laser intended for use.
 Use birefringent filter Part No. 0434-8974 for tuning near or slightly beyond 1100 nm. Align on the purple dot.

Table 8-3: Model 3900S Accessories and Replacement Parts

Description	Part Number	
Thin etalon assembly for <15 GHz operation	0445-0270	
Thick/thin etalon assembly for <1 GHz operation	0445-0050	
Model 3910 nitrogen purge box	Model 3910	
Model 3910 filter assembly	0449-7240	
Optional low %T output coupler	G0324-012 ¹	

To optimize performance at 980 nm with a 5 W pump source

Table 8-4: Hardware

Description	Part Number
Cooling assembly for Model 2060 and 2080 lasers ¹	0429-9200
Mirror spring (wavy washer) for M ₁ , M ₂ , and M ₃	0440-4430
Mirror spring (wavy washer) for M_4 , M_5 , M_6 , M_7 , and M_8	2801-4660
O-ring for M_4 , M_5 , M_6 , M_7 , and M_8	2504-2430
Threaded mirror assembly (holder, O-ring, wavy washer) for $\mathrm{M_4,M_5}$, and $\mathrm{M_6}$	0441-0160S
Threaded mirror assembly (holder, O-ring, wavy washer) for M_7 and M_8	0440-3600S
Nylon ferrules for ¼ in. tubing fitting	1300-0124
Polyethylene tubing, ¼ in., 50 ft.	2519-0120
Bulkhead fitting, ¼ in. tubing to ¼ in. tubing	1300-2014
Tool kit	0445-2460S
Routing mirror mounts, set of two with beam tubes ²	Model 341-04
Routing mirror, 488-514 nm broadband high reflector for use with Model 341-04	G0234-003

¹ Two included.² Optics not included.

Chapter 9

At Spectra-Physics, we take pride in the durability of our products. Considerable emphasis has been placed on controlled manufacturing methods and quality control throughout the manufacturing process. Nevertheless, even the finest precision instruments will need occasional service. We feel that our instruments have favorable service records compared to competitive products. We hope to demonstrate, in the long run, that we provide excellent service to our customers in two ways. First by providing the best equipment for the money, and second by offering service facilities that get your instrument back into action as soon as possible.

Spectra-Physics maintains major service centers in the United States, Europe, and Japan. Additionally, there are field service offices in major United States cities. When calling for service inside the United States, dial our toll-free number: 1 (800) 456-2552. To phone for service in other countries, refer to the Service Centers listing located at the end of this section.

Order replacement parts directly from Spectra-Physics. For ordering or shipping instructions, or for assistance of any kind, contact your nearest sales office or service center. You will need your instrument model and serial numbers available when you call. Service data or shipping instructions will be promptly supplied.

To order optional items or other system components, or for general sales assistance, dial 1 (800) SPL-LASER in the United States, or 1 (650) 961-2550 from anywhere else.

Warranty

All parts and assemblies manufactured by Spectra-Physics are unconditionally warranted to be free of defects in workmanship and materials for the period of time listed in the sales contract following delivery of the equipment to the F.O.B. point.

Liability under this warranty is limited to repairing, replacing, or giving credit for the purchase price of any equipment that proves defective during the warranty period, provided prior authorization for such return has been given by an authorized representative of Spectra-Physics. Warranty repairs or replacement equipment is warranted only for the remaining unexpired portion of the original warranty period applicable to the re paired or replaced equipment.

This warranty does not apply to any instrument or component not manufactured by Spectra-Physics. When products manufactured by others are

included in Spectra-Physics equipment, the original manufacturer's warranty is extended to Spectra-Physics customers. When products manufactured by others are used in conjunction with Spectra-Physics equipment, this warranty is extended only to the equipment manufactured by Spectra-Physics.

Spectra-Physics will provide at its expense all parts and labor and one way return shipping of the defective part or instrument (if required).

This warranty does not apply to equipment or components that, upon inspection by Spectra-Physics, discloses to be defective or unworkable due to abuse, mishandling, misuse, alteration, negligence, improper installation, unauthorized modification, damage in transit, or other causes beyond Spectra-Physics' control.

The above warranty is valid for units purchased and used in the United States only. Products with foreign destinations are subject to a warranty surcharge.

Warranty Return Procedure

Contact your nearest Spectra-Physics field sales office, service center, or local distributor for shipping instructions or an on-site service appointment. You are responsible for one-way shipment of the defective part or instrument to Spectra-Physics.

We encourage you to use the original Spectra-Physics packing boxes to secure instruments during shipment. If shipping boxes have been lost or destroyed, we recommend that you order new ones. Spectra-Physics will only return instruments in Spectra-Physics' containers.

Warning!



Always drain the cooling water from the laser head before shipping. Water expands as it freezes and will damage the laser. Even during warm spells or summer months, freezing may occur at high altitudes or in the cargo hold of aircraft. Such damage is excluded from warranty coverage.

Service Centers

Benelux

Telephone: (31) 40 265 99 59

France

Telephone: (33) 1-69 18 63 10

Germany and Export Countries*

Spectra-Physics GmbH Guerickeweg 7 D-64291 Darmstadt

Telephone: (49) 06151 708-0 Fax: (49) 06151 79102

Japan (East)

Spectra-Physics KK
East Regional Office
Daiwa-Nakameguro Building
4-6-1 Nakameguro
Meguro-ku, Tokyo 153

Telephone: (81) 3-3794-5511 Fax: (81) 3-3794-5510

Japan (West)

Spectra-Physics KK West Regional Office Nishi-honmachi Solar Building 3-1-43 Nishi-honmachi Nishi-ku, Osaka 550-0005

Telephone: (81) 6-4390-6770 Fax: (81) 6-4390-2760

e-mail: niwamuro@splasers.co.jp

United Kingdom

Telephone: (44) 1442-258100

United States and Export Countries**

Spectra-Physics

1330 Terra Bella Avenue Mountain View, CA 94043

Telephone: (800) 456-2552 (Service) or

(800) SPL-LASER (Sales) or (800) 775-5273 (Sales) or (650) 961-2550 (Operator)

Fax: (650) 964-3584

e-mail: service@splasers.com

sales@splasers.com

Internet: www.spectra-physics.com

^{*}And all European and Middle Eastern countries not included on this list.
***And all non-European or Middle Eastern countries not included on this list.

Appendix A

Material Safety Data Sheets (MSDS)

This section contains scanned copies of material safety data sheets (MSDS) that are supplied by the vendor and cover the various chemicals and compounds used in the Tsunami laser system, i.e., the compounds used in the filters in the *Model 3910* filter/purge unit.

Danger!



Read the MSDS sheets carefully before handling or disposing of filter elements. They may contain hazardous chemicals. Spectra-Physics has not independently determined the accuracy of the Material Safety Data Sheets, which were developed by the manufacturer of each chemical; therefore, we do not warrant the information contained therein. Dispose of these filter elements properly as indicated on the appropriate data sheet, and refer to your local environmental regulations regarding disposal. For further information, contact the chemical manufacturer at the address listed on each sheet.

Model 3900S CW, Ti:sapp	hire	Laser
-------------------------	------	-------

This page intentionally left blank.

Alltech MATERIAL SAFETY DATA SHEET

Alltech Associates, Inc. • 2051 Waukegan Road • Deerfield, IL 60015 • Phone: 708-948-8600 • Fax: 708-948-1078

STOCK NO:4034

PAGE: 1

ALSO APPLIES TO:*4034,4037,*4037

DATE:2/7/86

IDENTIFICATION

NAME: Borosilicate glass wool, silyated

REMARKS

The manufacturer of the glass wool used in these products indicates that an MSDS is not needed for their material. We have elected to provide some of the safety and common sense information we believe is useful.

The basic product has been treated with a silvation reagent to react with the silanol groups to make it more suitable for its intended application in chromatography. The glass wool is rendered hydrophobic. The excess reagent is removed and, to the best of knowledge, should not present any additional hazard.

CAS#:65997-17-3

SYNONYMS

Angel hair, fiberglass

TOXICITY HAZRDS

TOXICITY DATA

HEALTH HAZARD DATA

ACUTE EFFECTS

Overexposure: may cause temporary skin and upper respiratory irritation.

APPEARANCE AND ODOR

FIRST AID

Wash exposed areas with soap and warm water after handling.

FIRE AND EXPLOSION HAZARD DATA

EXTINGUISHING MEDIA

SPECIAL FIREFIGHTING PROCEDURES

UNUSUAL FIRE AND EXPLOSIONS HAZARDS

REACTIVITY DATA

INCOMPATIBILITIES STABILITY

Alltech MATERIAL SAFETY DATA SHEET (cont'd)

PAGE: 2

STOCK NO:4034 ALSO APPLIES TO:*4034,4037,*4037 DATE:2/7/86

HAZARDOUS COMBUSTION OR DECOMPOSITION PRODUCTS

HAZARDOUS POLYMERZATION WILL NOT OCCUR.

SPILL OR LEAK PROCEDURES

STEPS TO BE TAKEN IF MATERIAL IS RELEASED OR SPILLED

While wearing suitable protective devices, gather and sweep up. Place in plastic bags for disposal.

WASTE DISPOSAL METHOD

According to federal, state and local ordinances.

SPECIAL PROTECTION INFORMATION

SPECIAL PRECAUTIONS

Use respirator, such as 3M model 8710 or equivalent, for protection against nuisance dusts. Wear long sleeves and loose fitting clothing. Wash work clothes separartely from other clothing: then rinse washer thoroughly.

NOTICE: The information contained in the MSDS description is applicable exclusively to the chemical substances identified herein and for its intended use as an analytical reference standard or reagent and to the unit quantity intended for that purpose. The information does not relate to, and may not be appropriate for, any other applications or larger quantity of the substance described. Our products are intended for use by individuals possessing sufficient technical skill and qualification to use the material with suitable discretion and understanding of risk of handing any potentially hazardous chemical. The information has been obtained from sources believed to be reliable and accurate but has not been independently verified by Alltech Associates, Inc.

Accordingly, NO REPRESENTATION OR WARRANTY, EXPRESS OR IMPLIED, WITH RESPECT TO MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE is made with respect to the information contained herein.

ATTENTION: THIS PRODUCT IN TERMS OF CHEMICAL IDENTITY AND THE UNIT AMOUNT PROVIDED IS INTENDED FOR USE IN CHEMICAL ANALYSIS AND FOR HUMAN CONSUMPTION, NOR ANY OTHER PURPOSE.

Alltech MATERIAL SAFETY DATA SHEET

Alltech Associates, Inc. • 2051 Waukegan Road • Deerfield, IL 60015 • Phone: 708-948-8600 • Fax: 708-948-1078

STOCK NO:5779 PAGE:1 ALSO APPLIES TO:57720, 5779A, 8124, 5779 DATE:2/5/86

IDENTIFICATION

NAME: Activated Carbon

CAS#: 7740-44-0

SYNONYMS

Charcoal

TOXICITY HAZRDS

TOXICITY DATA

Activated carbon base material non-hazardous per Office of Materials Operation

HEALTH HAZARD DATA

ACUTE EFFECTS

Eyes: like any solid, mechanical irritation can occur.

APPEARANCE AND ODOR

Black irregular granules; free flowing when dry; without odor.

FIRST AID

Eyes: Flush with large volumes of water for 15 minutes. If irritation persists,

see a physician. Skin: Wash with soap and water; rinse well; dry.

FIRE AND EXPLOSION HAZARD DATA

EXTINGUISHING MEDIA

SPECIAL FIREFIGHTING PROCEDURES

UNUSUAL FIRE AND EXPLOSIONS HAZARDS

Dust can present a potential for dust explosion.

REACTIVITY DATA

INCOMPATIBILITIES

STABILITY

Stable.

AVOID

Close proximity to violate solvents; open flames; strong oxidizing or reducing

chemicals.

HAZARDOUS COMBUSTION OR DECOMPOSITION PRODUCTS

HAZARDOUS POLYMERZATION WILL NOT OCCUR.

SPILL OR LEAK PROCEDURES

Alltech MATERIAL SAFETY DATA SHEET (cont'd)

STOCK NO:5779 PAGE:2 ALSO APPLIES TO:57720, 5779A, 8124, 5779 DATE:2/5/86

STEPS TO BE TAKEN IF MATERIAL IS RELEASED OR SPILLED

Collect with broom and shovel. Wear a particle mask, goggles and gloves. If media becomes wet, store for waste disposal in plastic-lined containers. Do not enter large closed vessels without a self-contained breathing apparatus-especially when carbon is wet as wet activated carbon removes oxygen from air.

WASTE DISPOSAL METHOD

Dispose of unspent carbon in refuse container.

SPECIAL PROTECTION INFORMATION

Respiratory: Dust or particle mask.

Mechanical: Dust filter

Protective Gloves:Corrosive chemical resistant

SPECIAL PRECAUTIONS

While handling, use a particle mask, goggles and gloves. Store away from violate organic solvents and moisture (to preserve media)

NOTICE: The information contained in the MSDS description is applicable exclusively to the chemical substances identified herein and for its intended use as an analytical reference standard or reagent and to the unit quantity intended for that purpose. The information does not relate to, and may not be appropriate for, any other applications or larger quantity of the substance described. Our products are intended for use by individuals possessing sufficient technical skill and qualification to use the material with suitable discretion and understanding of risk of handing any potentially hazardous chemical. The information has been obtained from sources believed to be reliable and accurate but has not been independently verified by Alltech Associates, Inc.

Accordingly, NO REPRESENTATION OR WARRANTY, EXPRESS OR IMPLIED, WITH RESPECT TO MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE is made with respect to the information contained herein.

ATTENTION: THIS PRODUCT IN TERMS OF CHEMICAL IDENTITY AND THE UNIT AMOUNT PROVIDED IS INTENDED FOR USE IN CHEMICAL ANALYSIS AND FOR HUMAN CONSUMPTION, NOR ANY OTHER PURPOSE.



MATERIAL SAFETY DATA SHEET

(Essentially similar to U. S. Department of Labor Form OSHA-20)

Do Not Duplicate This Form. Request an Original

	31 (9 - 2) 31 (9 - 2) 31 (9 - 2) 31 (9 - 2)	I. PRODUCT.II	ENTIFICA	TION			
PRODUCT Molecular Siev	e Type 5A					·	
CHEMICAL NAME Mo	lecular Siev	e Type 5A					
SYNONYMS Zeolites			CHEMICAL FAMILY Calcium Sodium Alumino Silicates				
FORMULA CaNa SiO ₂ AL ₂ O ₃			MOLECULAR WEIGHT				
TRADE NAME AND SYNOR	NYMS	Linde Molecular S	eve and MOL	SIV			
TO SERVICE OF THE SER	3.4.21(4) 4",	II HAZARDOU	SINGRED	IENTS			
No TLV's have b		hed for this product			000 (1976)	or ACGIH	1977
	T	ACGIH (1977)		MATERIAL		Wt (%)	ACGIH (1977)
MATERIAL	Wt (%)	TLV-TWA (Units)		MATERIAL		VV (76)	TLV-TWA (Units)
l							
			ICAL DAT				
		III. PHYS	AN HOLL	FREEZIN	C POINT	144	
BOILING POINT, 760 mm.		NA			RESSURE	A T 20°C	NA
SPECIFIC GRAVITY (H2O	= 1)	1.1		SOLUBIL		Negli	· · · · · · · · · · · · · · · · · · ·
VAPOR DENSITY (air = 1)	NA			WATER,			
PER CENT VOLATILES BY VOLUME	NA			(Butyl Ac	etate = 1)		NA ————————————————————————————————————
APPEARANCE AND ODOR		dor — depending on p					
	i, IV.	FIRE AND EXPL	OSION HA	ZARD DA	TA:		
	onflammabl		AU1	TOIGNITION PERATURE	N P		
FLAMMABLE LIMITS IN AIR, % by volume		LOWER	NA		UPPER	NA	
EXTINGUISHING MEDIA	١	NA					
SPECIAL FIRE FIGHTING PROCEDURES	١	NA					
UNUSUAL FIRE AND EXPLOSION HAZARDS	1	None known					
	\$ 2 to \$ 10 to	EMERGENCY	PHONE NU	MBER		$E^{(r)}$	
IN CASE OF EMERGE	NCIES invol	lying this material, furth		is available at a	all times at th	is telephone	number:
	For	304: routine information co	744-3487 ntact your local	Linde Supplie	r.		
While Union Carbide Corporation regarding the tests conducted, the responsibility. They are offered							
responsibility. They are offered determined by the user to be in ac	cordance wit	h applicable Federal, St	ate, and local la	ws and regulat	lons.		

UNION CARBIDE CORPORATION . LINDE DIVISION . 270 PARK AVENUE, NEW YORK, N.Y. 10017


Molecular Sieve Type 5A PRODUCT: V-HEALTH HAZARD PATAS See Section II THRESHOLD LIMIT VALUE EFFECTS OF OVEREXPOSURE AND EMERGENCY AND FIRST AID PROCEDURES Irritation of eyes, nose and throat by dust. In case of eye contact, immediately flush eyes with plenty of water for at least 15 minutes. If irritation persists, see a physician. STABILITY UNSTABLE STABLE CONDITIONS TO AVOID None known Sudden contact with high concentrations of chemicals having high heats of adsorption INCOMPATIBILITY (materials to avoid) such as Olefins, HCL, etc. HAZARDOUS DECOMPOSITION PRODUCTS None known HAZARDOUS POLYMERIZATION CONDITIONS TO AVOID May Occur | Will not Occur None known VII. SPILL OR LEAK PROCEDURES STEPS TO BE TAKEN IF MATERIAL IS RELEASED OR SPILLED Sweep up and place in a waste disposal container. Flush area with water. Avoid raising dust. WASTE DISPOSAL METHOD Bury in a landfill. VIII. SPECIAL PROTECTION INFORMATION RESPIRATORY PROTECTION (specify type) If there is excessive dustiness, wear a respirator selected as per OSHA 29 CFR 1910.134 SPECIAL LOCAL EXHAUST See OSHA 29 CFR 1910.134 None **VENTILATION** OTHER MECHANICAL (general) None See OSHA 29 CFR 1910.134 EYE PROTECTION - Safety glasses PROTECTIVE GLOVES Not required but recommended or goggles selected as per OSHA 29 CFR 1910.133 OTHER PROTECTIVE EQUIPMENT IX. SPECIAL PRECAUTIONS Causes eye irritation. Breathing dust may be harmful. May cause skin irritation. Open container slowly to avoid dust. Do not get in eyes. Avoid breathing dust and prolonged contact with skin. Use with adequate ventilation. Keep container closed. Wash thoroughly after handling. Do not ingest. OTHER HANDLING AND STORAGE CONDITIONS


F-4503 80-0946 11/80 5C


pH Range if in Aqueous Slurry - 8 - 11


Lithographed in U.S.A.


Notes














Report Form for Problems and Solutions

We have provided this form to encourage you to tell us about any difficulties you have experienced in using your Spectra-Physics instrument or its manual—problems that did not require a formal call or letter to our service department, but that you feel should be remedied. We are always interested in improving our products and manuals, and we appreciate all suggestions. Thank you.

rioiii.	
Name	
Company or Institution	
Department	
Address	
Instrument Model Number	
Instrument Model Number Problem:	
Suggested Solution(s):	

Mail To:

Erom:

Spectra-Physics, Inc. ISL Quality Manager 1330 Terra Bella Avenue, M/S 15-50 Post Office Box 7013 Mountain View, CA 94039-7013 U.S.A.

E-mail: sales@splasers.com www.spectra-physics.com

FAX to:

Attention: ISL Quality Manager (650) 961-7101